Coaxial Sleeve

Gebruuchsanweisung

Besondere Beachtung

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Instructions for Use

Special notice

 Please read the following information and operating instructions carefully.

 **Caution:** Federal law restricts this device to sale by or on the order of a physician.

The device may only be used by qualified medical staff in accordance with these user instructions.

PAJUNK® does not recommend any particular treatment method. Professional medical staff are responsible for the way in which the device is used and for patient selection.

Failure to comply with the user instructions invalidates the warranty and puts patient safety at risk.

If used in combination with other products, it is essential that the compatibility information and user instructions for these other products are taken into account. A decision regarding the combined use of devices from different manufacturers (where they do not constitute treatment units) is the responsibility of the user.

 The device must not be used under any circumstances if there are good reasons to suspect incompleteness, damage or loss of sterility.

 Only devices in perfect condition, which are within the sterile expiry date marked on the label, in undamaged packaging, may be used.

Device description/ compatibility


 Please see the current declaration of conformity for product numbers and the scope of these instructions for use.

Accessory for PAJUNK® instruments for soft tissue biopsy, to be used solely with the PAJUNK® biopsy systems DeltaCut, CoreCut and PrimoCut.

 In addition to these instructions for use, the relevant information also applies according to the corresponding specialist literature and current state of the art and knowledge.

 Users must provide patients with information.

 **Warning:** The biopsy cannula is not suitable for MRI use!

Intended use

Accessory for the removal of tissue samples from soft tissue.

Indikation

Cytological and/ or histological examination.

Contraindications

No therapeutic consequence, uncooperative patient, ascites, poorly visible organs, severe coagulation disorder, unsafe access route, aneurysm, phaeochromocytoma, Echinococcus, injury to neighbouring organs (lung injury, gall-bladder injury, bowel injury), infections, hypersensitivity reactions to the local anesthetic, cardiovascular disorders following administration of analgesics or sedatives.

Under no circumstances is the device to be used in the event of known material incompatibilities and/ or known interactions.

Complications

Incorrect puncture, coagulation disorder, poor general health, haematoma in the region of the target area, pneumothorax, haemothorax, vascular injury, arterio-biliary fistula.

If complications occur while using the device, follow the protocols of your organisation. If this does not resolve the complications, or if they are regarded as serious or untreatable, carefully stop the procedure and remove invasive device components from the patient.

Warnings

for sterile product:

This is a disposable medical device for use with only one patient.

This device must not be re-used under any circumstances.

This device must not be resterilised under any circumstances.

The materials used in the manufacture of this device are not suitable for reprocessing or resterilisation.

This device is not designed to be reprocessed or resterilised.

-  **Unauthorised re-use or reprocessing**
- ... can cause the device to lose the performance properties intended by the manufacturer.
 - ... leads to a significant risk of cross-infection/ contamination as a result of potentially inadequate processing methods.
 - ... may cause the device to lose significant functional properties.
 - ... may cause materials to break down and lead to endotoxic reactions caused by the residues.

 **in the application:**

1. When using the biopsy system make sure that the biopsy cannula does not become bent.
2. It is vital to comply with good clinical practice. Serious post-procedural complications include deep wound infections requiring extensive surgery.
3. The biopsy specimen must only be taken in a clinical environment.
4. Ensure that you take appropriate measures before the puncture to safeguard the biopsy specimen for haematopathological examination.

 **regarding puncture:**

1. If you unexpectedly come into contact with bone, change the direction of the cannula. Do not try to overcome bone resistance. Failure to adhere to these rules could cause the cannula to bend or break.
2. Take care to use devices of suitable dimensions (diameter, length), especially when treating obese patients and children.

 **for injection:**

Always ensure that the injection site is aseptic.

 **for use with other compatible products:**

When using multiple components, familiarise yourself with their operation before use by checking connections and passages (cannulas, adapters).

Further warning indications

1.  **Caution: Sharp object warning.** The device or device components may, depending on the type of tip, have sharp edges or tips. Various infectious pathogens can be transmitted if a stab wound occurs. For practical purposes, the most important of these are the human immunodeficiency virus (HIV), the hepatitis B virus (HBV) and the hepatitis C virus (HCV).
2. You must routinely take general precautions for handling blood and bodily fluids when using and disposing of the device, due to the risk of contact with blood-borne pathogens.

Sequence of use

 When performing a biopsy, follow the correct procedure as laid down in the instructions for use of the biopsy system you are using.

1. Use suitable, sterile techniques when performing the biopsy.
2. Before positioning the system, make a stab incision in the skin, where necessary, to facilitate penetration.
3. Pass the coaxial sleeve through the incision, with the help of an imaging system if necessary, until the tip of the coaxial sleeve is proximal to the area from which the biopsy is to be taken.
4. A DeltaCut, CoreCut or PrimoCut biopsy cannula can now be inserted. The procedure is described in the relevant instructions for use.

Proceed as follows after obtaining the biopsy specimen:

5. Remove the instrument and/ or the cannula from the coaxial sleeve.
6. If further biopsy samples are to be taken, the coaxial sleeve can be left at the puncture site. It is vital that the conditions around the puncture site are hygienic and aseptic.
7. Once the procedure has been completed, carefully remove the coaxial sleeve.
8. Treat the incision site.

Use and storage conditions

	Temperature limit	+10 °C to +30 °C
	Humidity limitation	20 % to 65 %
	Keep away from sunlight	
	Keep dry	

General information

The devices are manufactured in accordance with globally applicable guidelines for hazardous substances.

 **PAJUNK®** GmbH Medizintechnologie, Karl-Hall-Strasse 1, 78187 Geisingen, Germany.

Key to symbols used in labelling

	Manufacturer		Consult instructions for use
	Use-by date		Caution: Federal law restricts this device to sale by or on the order of a physician
	Catalogue number		MR unsafe
	Sterilized using ethylene oxide		Advice
	Do not re-sterilize		Information
	Do not use if package is damaged		Product is in conformity with the applicable requirements set out in Community harmonization legislation and is monitored by a notified body
	Keep dry		Sharp object warning
	Humidity limitation		Does not contain Phthalates (acc. to sec. 7.5 of Annex I 93/42/EWG)
	Do not re-use		Natural rubber has not been used as a component in the manufacture of this product
	Caution		Quantity
	Date of manufacture		
	Batch code		
	Keep away from sunlight		
	Temperature limit		

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PAJUNK® GmbH Medizintechnologie, Karl-Hall-Strasse 1, 78187 Geisingen, Alemania.
PAJUNK[®] GmbH