SEE DETAILS

THAT MAKE

bk3000

THE DIFFERENCE

- Dedicated urology ultrasound system with a simple workflow and small footprint
- Wide range of sterilizable urology transducers

bk)

medical

- Advanced modes, including Elastography and Contrast Imaging, to help identify regions of interest
- Immediate, auto-optimized imaging using No-touch Autogain and Auto Focus features

as specialized as you

UNPARELLELED IMAGING TO GUIDE YOUR PROCEDURES

High-resolution images with superb anatomical details are enabled by faster processing, advanced auto-optimization algorithms, and BK's high performance transducers.

TRANSRECTAL PROSTATE BIOPSY





lesions and perform targeted biopsies.

FUSION



procedure.

Reduce registration errors and procedure steps with Predictive Fusion.





UROGYNECOLOGY



KIDNEY/BLADDER





(9085) 6C2 (9040) rved Curved 9C2 (9002) 6C2s (9023) Curved Small Curved

Right kidney resistive index measurement 6C2s transducer

Curved array transducers provide optimal image quality for advanced abdominal kidney imaging, measurement, and biopsy.

PENILE/TESTES



DEDICATED UROLOGY ULTRASOUND SOLUTION

DESIGNED FOR UROLOGY

- A small footprint for a variety of clinical settings.
- Sensitive color Doppler provides excellent real-time visualization of vascular perfusion.
- Advanced imaging modes, including elastography and contrast imaging, help identify regions of interest.
- Rapid boot time and up to two hours of plug-free imaging.
- Sterilizable remote control for control of the system in the sterile field.

FULLY STERILIZABLE, DEDICATED UROLOGY TRANSDUCERS

- Ergonimically designed for a secure and comfortable grip.
- Long transducer cables improve flexibility in the OR or clinic.
- Most transducers on the bk3000/bk5000 feature a Smart[™] button that enables you to activate the transducer, freeze, store or print images at the press of a button.



A selection of BK Medical's specialized transducers

SINGLE-USE NEEDLE GUIDES

• Single-use needle guides help you perform reliable puncture and biopsy procedures quickly.



Simultaneous Biplane single-use needle guide



Endfire single-use needle guide



A selection of BK Medical single-use needle guides

BK Medical

8 Centennial Drive Peabody MA 01960 USA T +1 978 326 1300 bkmedical.com North America Sales & Service BK Medical 8 Centennial Drive Peabody, MA 01960 USA T +1 978 326 1300 F +1 978 326 1399 bkmedical.com

bk3000/bk5000

Europe and Rest of World Sales, Service & Design Center BK Medical Mileparken 34 2730 Herlev Denmark T +45 4452 8100 F +45 4452 8199 bkmedical.com

