

# THE IPAS MANUAL VACUUM ASPIRATION TECHNOLOGY PRODUCT LINE CATALOGUE

A COST-EFFECTIVE SOLUTION FOR UTERINE EVACUATION<sup>1</sup>

- BACKED BY 40 YEARS OF GLOBAL MARKET USE AND IPAS' LONGSTANDING QUALITY REPUTATION
- FDA-LISTED AND DESIGNED AND MANUFACTURED CONFORMING TO ISO 13485
- SMALL, INEXPENSIVE, PORTABLE, QUIET, AND SAFE<sup>2,3,4</sup>
- IDEAL FOR PERFORMING PROCEDURES IN OUTPATIENT SETTING
- IDEAL FOR SETTINGS WITH NO ELECTRICITY OR WITH SPACE CONSTRAINTS
- LIMITED INVESTMENT REQUIRED
- EASY TO USE

## CLINICAL BENEFITS IN MANUAL UTERINE EVACUATION

	IPAS ASPIRATORS	EVA	SHARP CURETTAGE
SAFE	X	X	
LOWER INCIDENCE OF SIDE EFFECTS/COMPLICATIONS <sup>5</sup>	X	X	
COST EFFECTIVE	X		X*
PORTABLE	X		X
QUIET	X		X
SINGLE USE	X		
REUSABLE	X	X	X
OUTPATIENT/NO ANESTHESIA <sup>6</sup>	X	X	
QUICK RECOVERY TIME	X	X	
NO ELECTRICITY REQUIRED	X		X
REQUIRES LIMITED SPACE	X		X
REQUIRES LIMITED INVESTMENT	X		X

\*When used without anesthesia

#### References:

1. Veruyll DAA, Crowther CA. (1993). Suction v. conventional curettage in complete abortion: a randomized clinical trial. S Afr Med J.83:13-15.
2. Blumenthal, Paul D and R.E. Remsburg. 1994. A time and cost analysis of the management of incomplete abortion with manual vacuum aspiration. Int J Gynecol Obstet, 1994, 45: 261-267.
3. De Jonge ETM, Pattinson, RC, Makin JD et al. (1994). Is ward evacuation for uncomplicated incomplete abortion under systematic analgesia safe and effective: a randomized clinical trial. S Afr Med J, 84:481-483.
4. Fawcus S, McIntyre J, Jewkes RK et al. (1997). Management of incomplete abortions at South African public hospitals. S Afr Med J, 87:438-442.
5. World Health Organization, Department of Reproductive Health and Research. Safe Abortion: Technical and policy guidance for health systems. Second edition. 2012. Annex 5, Recommendation 1, page 113.
6. Hale RW, Kobara TY, Sharma SD, Tsuei JJ, Gramlich EP, Nakayama RT. Office termination

WomanCare Global Trading CIC is a UK Registered Community Interest Company (No.07231920).  
E-mail: [contact@dktwomancare.org](mailto:contact@dktwomancare.org)  
[www.dktwomancare.org](http://www.dktwomancare.org)

## PRODUCT PORTFOLIO

### ASPIRATORS

- Ipas MVA Plus® Aspirator
- Ipas Double Valve Aspirator
- Ipas Single Valve Aspirator

### CANNULAE

- Ipas EasyGrip® Cannulae
- Flexible Karman Cannulae

### DILATORS AND ACCESSORIES

- Adapter sets, dilators and aspirator accessories also available

### KITS

- Miscarriage Management Kit

### TRAINING TOOLS

- Pelvic Model

WomanCare GLOBAL A COMPANY OF **dkt** INTERNATIONAL

An ISO 9001-certified company

# Ipas Double-Valve Aspirator

For uterine evacuation and endometrial biopsy

EFFECTIVE · DISPOSABLE · AFFORDABLE · VERSATILE



## CLINICAL INDICATIONS

- Induced or incomplete abortion up to 12 weeks since the woman's last menstrual period
- Miscarriage management
- Endometrial biopsy

## PRODUCT CHARACTERISTICS

- Individually packed
- Provides a holding capacity of 60 cc
- Provides a minimum vacuum of 24 inches
- Maintains vacuum for a minimum of 30 minutes
- Cylinder marks in 10 cc increments to measure evacuated uterine contents
- Compatible with Ipas EasyGrip® Cannulae and the 12 mm Flexible Karman Cannula
- Latex-free plastic
- Shelf life: 5 years from month of manufacture

## ALSO AVAILABLE FROM WOMANCARE GLOBAL



### IPAS MVA PLUS® ASPIRATOR

For uterine evacuation  
Product Code: PLUS



### IPAS EASYGRIP® CANNULAE

Product Code: SR4 to SR10, and SR12  
or SR4-SU to SR10-SU, and SR12-SU



### FLEXIBLE KARMAN CANNULAE

For uterine evacuation  
Product Code: FC4 to FC10, and FC12