

# THE IPAS MANUAL VACUUM ASPIRATION TECHNOLOGY PRODUCT LINE CATALOGUE

A COST-EFFECTIVE SOLUTION FOR UTERINE EVACUATION<sup>1</sup>

- BACKED BY 40 YEARS OF GLOBAL MARKET USE AND IPAS' LONGSTANDING QUALITY REPUTATION
- FDA-LISTED AND DESIGNED AND MANUFACTURED CONFORMING TO ISO 13485
- SMALL, INEXPENSIVE, PORTABLE, QUIET, AND SAFE<sup>2,3,4</sup>
- IDEAL FOR PERFORMING PROCEDURES IN OUTPATIENT SETTING
- IDEAL FOR SETTINGS WITH NO ELECTRICITY OR WITH SPACE CONSTRAINTS
- LIMITED INVESTMENT REQUIRED
- EASY TO USE

## CLINICAL BENEFITS IN MANUAL UTERINE EVACUATION

	IPAS ASPIRATORS	EVA	SHARP CURETTAGE
SAFE	X	X	
LOWER INCIDENCE OF SIDE EFFECTS/COMPLICATIONS <sup>5</sup>	X	X	
COST EFFECTIVE	X		X*
PORTABLE	X		X
QUIET	X		X
SINGLE USE	X		
REUSABLE	X	X	X
OUTPATIENT/NO ANESTHESIA <sup>6</sup>	X	X	
QUICK RECOVERY TIME	X	X	
NO ELECTRICITY REQUIRED	X		X
REQUIRES LIMITED SPACE	X		X
REQUIRES LIMITED INVESTMENT	X		X

\*When used without anesthesia

#### References:

1. Veruyil DAA, Crowther CA. (1993). Suction v. conventional curettage in complete abortion: a randomized clinical trial. S Afr Med J.83:13-15.
2. Blumenthal, Paul D and R.E. Remsburg. 1994. A time and cost analysis of the management of incomplete abortion with manual vacuum aspiration. Int J Gynecol Obstet, 1994, 45: 261-267.
3. De Jonge ETM, Pattinson, RC, Makin JD et al. (1994). Is ward evacuation for uncomplicated incomplete abortion under systematic analgesia safe and effective: a randomized clinical trial. S Afr Med J, 84:481-483.
4. Fawcus S, McIntyre J, Jewkes RK et al. (1997). Management of incomplete abortions at South African public hospitals. S Afr Med J, 87:438-442.
5. World Health Organization, Department of Reproductive Health and Research. Safe Abortion: Technical and policy guidance for health systems. Second edition. 2012. Annex 5, Recommendation 1, page 113.
6. Hale RW, Kobara TY, Sharma SD, Tsuei JJ, Gramlich EP, Nakayama RT. Office termination

WomanCare Global Trading CIC is a UK Registered Community Interest Company (No.07231920).  
E-mail: [contact@dktwomancare.org](mailto:contact@dktwomancare.org)  
[www.dktwomancare.org](http://www.dktwomancare.org)

## PRODUCT PORTFOLIO

### ASPIRATORS

- Ipas MVA Plus® Aspirator
- Ipas Double Valve Aspirator
- Ipas Single Valve Aspirator

### CANNULAE

- Ipas EasyGrip® Cannulae
- Flexible Karman Cannulae

### DILATORS AND ACCESSORIES

- Adapter sets, dilators and aspirator accessories also available

### KITS

- Miscarriage Management Kit

### TRAINING TOOLS

- Pelvic Model


WomanCare GLOBAL A COMPANY OF **dkt** INTERNATIONAL

An ISO 9001-certified company

# IPAS EasyGrip® Cannulae

The tactile response of a rigid curette with the gentle probe of a flexible cannulae

SEMI-RIGID · INNOVATIVE · VERSATILE · REUSABLE



## PRODUCT CHARACTERISTICS

- Ready to use and sterile packed
- Reusable up to 25 times
- Permanently-integrated winged base
- Color-coded by size
- Uterine depth markers from 6 to 11 cm printed on cannulae
- Sized 4 mm to 8 mm - two "whistle cut" openings
- Sized 9 mm, 10 mm, and 12 mm - one large "scoop cut" opening
- Approximately 24 cm (9.5 in) in length
- Latex-free plastic
- Shelf life: 5 years from month of manufacture

### COLOR CODE



### PRODUCT CODE/DIAMETER CHART

Size	4 mm	5 mm	6 mm	7 mm	8 mm	9 mm	10 mm	12 mm
Code	SR4	SR5	SR6	SR7	SR8	SR9	SR10	SR12

## ALSO AVAILABLE FROM WOMANCARE GLOBAL



### IPAS DOUBLE VALVE ASPIRATOR

For uterine evacuation  
Product Code: DVS-10 or DVS-SU



### IPAS MVA PLUS® ASPIRATOR

For uterine evacuation  
Product Code: PLUS



### UTERINE PELVIC MODEL

For uterine evacuation training  
Product Code: UPM