


Salvafast®

Restraint and positioning
Contención y posicionamiento



QUICK
RETENTION SYSTEM

SISTEMA DE RETENCIÓN
RÁPIDA

Salvafast Transfer · Salvafast

Salvafast

TRANSFER

For
STRETCHER
Para
CAMILLA

Patent ES1133355U

QUICK RETENTION AND TRANSFER SYSTEM FOR STRETCHER

Innovative system for fast and safe restraint and transfer of disturbed, violent, aggressive, or disoriented patients for emergency stretchers. It allows to transfer the patient to another stretcher or bed without releasing the restraint. It consists of a sturdy sliding base that includes a set of harness, wrist, and ankle straps with a quick release double Velcro® fastening system.



SISTEMA DE RETENCIÓN RÁPIDA Y TRANSFERENCIA PARA CAMILLA

Innovador sistema para una rápida y segura retención y transferencia de personas alteradas, violentas, agresivas o desorientadas, para camillas de emergencia. Permite, en situaciones de emergencia, la transferencia de la persona a otra camilla o cama sin necesidad de liberar la contención. Está formado por una resistente base deslizante con un conjunto de arnés, muñequeras y tobilleras con sistema rápido de cierre de doble Velcro®.

	Measures/Medidas	Code/Código	Size/Talla
Base	172 x 45 cm	SVT3710	One/Única
Wristbands/Muñequeras	(P) 15 - 30 cm		
Ankle bands/Tobilleras	(P) 20 - 40 cm		



EN 1865
INTERNATIONAL
COMPLIANT



MONOBLOCK



- Harness + wristbands + ankle bands in one piece
- Designed specifically for stretchers
- Improved versatility: quick restraint + transfer
- Colored bands to identify the function of each one:

BLUE: fixation to the stretcher

BLACK: patient retention

ORANGE: handles for transfer

- Quick and safe application

- Easy to sanitize and store

- Recommended for patients up to 150 kg

- Arnés + muñequeras + tobilleras en una única pieza

- Diseño específico para camillas

- Mayor versatilidad: contención rápida + transferencia

- Cintas de color para identificar la función de cada una:

AZULES: fijación a la camilla

NEGRAS: retención del paciente

NARANJAS: asas para transferencia

- Aplicación rápida y segura

- Fácil de higienizar y almacenar

- Recomendado para pacientes hasta 150 kg

Harness

Arnés



Belt

Cinturón



Wristbands

Muñequeras



Ankle bands

Tobilleras



Handles for transfer

Asas para transferencia



Sliding base

Base deslizante



Comfortable and soft touch

Confortable y tacto suave



Compact storage bag

Bolsa para almacenamiento compacto

